

Honey - Healing without the Sting!

SanoSkin® Melladerm Honey Wound Gel

Sanoskin® Melladerm® honey wound gel is specially formulated and effective for the treatment of oncologic and/or acute and surgical wounds (contaminated) chronic wounds deep and superficial wounds with low to high exudates leg ulcers and pressure sores burns

Key Benefits

- Offers a more gentle approach than pure honey
- Maintains moist wound healing
- Utilises the antibacterial strengths of honey

Key Properties

When the gel is in contact with the wound, fluid is extracted from the surrounding tissues through osmosis creating a moist healing environment. Moist wound healing is known to optimise the wound healing process because this is an ideal physiological environment for the cells⁽¹⁾. This environment also

- **facilitates autolysis**
- **promotes better epithelial cell migration**
- **is an excellent debriding agent (autolysis)**
- **is painless in contrast to other honey products**
- **reduces malodour**
- **will dilute gradually and can be removed easily with a wound cleanser**



Sanoskin® Melladerm® is a primary wound dressing that can be combined with most commonly used secondary dressings, such as Sanoskin® foam or Sanoskin® NET plus many other wound dressings and will fit easily into any wound protocol.

Sanoskin® Melladerm® can be applied following light therapy treatment with the Sanoskin® Photizo® LED Light Therapy Device.

Directions for Use

For single patient use only.

Remove lid, clean the mouth of the tube and the top of the lid with an alcohol swab.

Push the top of the lid into the silver coating to break the seal.

Squeeze the tube and fill the wound, taking care not to touch the wound.

Use a (sterile) spatula when necessary.

Place a dressing on top of the gel to cover the wound.

Sanoskin® Melladerm® Plus can stay on the wound for up to 48 hours on average.

Renew the gel in the wound as soon as the gel has disappeared.

Store at room temperature, maximum 25°C / 77°F.

Order Code	Size	PIP Code
MD50	50g tube	334-7762

Contraindications No contraindications are known to date. Do not use on individuals with a known sensitivity to the gel or its components. Do not use if package is damaged or opened. Do not use the gel in cases where the patient is very susceptible to an infection e.g. in very large burn wounds.

1. Vanessa Jones, Keith Harding. Moist Wound Healing In: Krasner DL, Roseheaver GT, Sibbald RG, eds Chronic Wound Care: A Clinical Source Book for Healthcare Professionals. 3rd ed. Wayne, Pa: HMP Communications, 2001:245-252